


## Intrauterine aspiration





## INTRACURETTE SR

Curetta semirigida tipo "Karman"  
per aspirazione intrauterina

Flexible intrauterine suction curette - "Karman" style

Cod.	Descrizione - Description	Colore Colour
720156	lung. 24 cm - length 24 cm Ø esterno 4 mm - O.D. 4 mm	bianco white
720157	lung. 24 cm - length 24 cm Ø esterno 5 mm - O.D. 5 mm	beige beige
720114	lung. 24 cm - length 24 cm Ø esterno 6 mm - O.D. 6 mm	grigio grey
720115	lung. 24 cm - length 24 cm Ø esterno 7 mm - O.D. 7 mm	giallo yellow
720116	lung. 24 cm - length 24 cm Ø esterno 8 mm - O.D. 8 mm	arancio orange
720117	lung. 24 cm - length 24 cm Ø esterno 9 mm - O.D. 9 mm	rosso red
720118	lung. 24 cm - length 24 cm Ø esterno 10 mm - O.D. 10 mm	verde green
720119	lung. 24 cm - length 24 cm Ø esterno 11 mm - O.D. 11 mm	azzurro cyan
720120	lung. 24 cm - length 24 cm Ø esterno 12 mm - O.D. 12 mm	blu blue

 Per interruzioni volontarie di gravidanza o a seguito di aborto spontaneo. La flessibilità del materiale permette un curettaggio delicato dopo l'evacuazione. Disponibile in 9 diverse misure con codice colore per una facile identificazione. Si connette direttamente al tubo di aspirazione.

 For uterine evacuation and curettage in voluntary pregnancy termination or following miscarriage. Flexibility of the material allows a gentle curettage after evacuation. It is available in 9 different colour-coded sizes for an easy identification. It can be connected directly to the suction tubing.

Scatola da 50 pezzi.  
Fornito in confezione singola, monouso, sterilizzato a EtO e marcato CE.

50 pcs. per box.  
Supplied in single pouch, single use, EtO sterilized and CE marked.

