




**Sterilcontainer
Sterilizing Containers
Contenedores de Esterilización
Container di Sterilizzazione
Conteneurs de Stérilisation**





Benefits:

- Long lasting sterilization.
- Maximum flexibility in the preparation of surgical sets and sterilizable material.
- Protects surgical instruments from damage.
- Standard size allows easy stacking of containers
- Provides safe transportation of surgical instruments.
- Risk free return to storage against contamination.
- Cost-effective.

Features:

All dimensions comply with DIN 58952 standard, thus guaranteeing an economical rate of utilization of the sterilizers.

Container lids and bottoms are made of anodized aluminum-alloy, therefore, due to the thermal conductivity of the material an optimum drying phase is achieved.

All standard size containers can be stacked, thus space and economic storage problems are eliminated.

Wire baskets can also be stacked to prevent valuable instruments from being damaged during preparation of the sets.

Removal containers; ideal for safe return of dirty instruments and textiles without any risk of contamination during transportation after operation.

Wrapping drapes provide a safe sterile medium in the containers and wire baskets; and guarantee the sterilized contents against contamination when they are taken out for operation.

Ensures an economic and long lasting sterilization.

Sterilization:

We recommend following steam sterilization procedures for our containers with perforated lids and non-perforated bottoms, as specified in DIN 58946/6 standard.

For fractional flow procedure:

A few jets of saturated steam of $p_e \geq 1$ bar is generally sufficient to remove the air out of sterilizing chamber.

For fractional vacuum procedure:


















- Repeated evacuation up to a pressure of < 130 mbar in exchange with steam inflow up to a pressure above or below an atmospheric pressure.
- Maintaining of steam inflow until the working overpressure is reached.

For steam injection procedure (Vacuum procedure):

- Single evacuation up to a pressure of < 70 mbar with simultaneous low steam inflow.
- Maintaining of steam inflow until the working pressure is reached.

For further assistance regarding all other gravitation and prevacuum procedures as specified in DIN 58953/9 standard, please refer to the table on next page.

Suitability of our sterile containers for steam sterilization procedures ¹⁾ as per DIN 58953 part 9 (issue October 2002)

Sterilization container	Packing		Steam sterilization procedure acc. to DIN 58946, part I				
	Description	Symbolic representation and packing system	Gravitation procedure	Pre-vacuum procedure	Fractional vacuum procedure	Fractional flow procedure	Steam injection procedure
with filters (StB-F)	Textiles in sterilization container with filter in lid and bottom acc. to DIN 58952 part I		+	+	+	+	+
			○	○	+	+	+
	Textiles in sterilization container with filter in lid, bottom not perforated acc. to DIN 58952 part I		-	-	+	+	+
			-	-	+	+	○
	Textiles in sterilization container with filter in bottom, lid not perforated acc. to DIN 58952 part I		○	+	+	+	+
			-	○	+	+	+
	Instruments in wire basket acc. to DIN 58952 part III in sterilization container with filter in lid and bottom acc. to DIN 58952 part I		+ ³⁾	+	+	+	+
	Instruments in wire basket acc. to DIN 58952 part III, in sterilization container with filter in lid, bottom not perforated acc. to DIN 58952 part I		-	+ ²⁾	+ ²⁾	+ ²⁾	+ ²⁾
	Instruments in wire basket acc. to DIN 58952 part III, in sterilization container with filter in bottom acc. to DIN 58952 part I		○	+	+	+	+
			-	-	+	+	+
with Valves (StB-V)	Textiles in sterilization container with valves in lid and bottom acc. to DIN 58952 part I		-	○	+	+	○
			-	-	+	+	○
	Textiles in sterilization container with valves in lid no valves on the bottom		-	-	+	+	-
			-	-	+	+	-
	Instruments in wire basket in sterilization container with valve in lid and bottom		-	+ ²⁾	+ ²⁾	+ ²⁾	+ ²⁾
	Instruments in wire basket in sterilization container with valve in lid not perforated		-	+ ²⁾	+ ²⁾	+ ²⁾	+ ²⁾
			-	+ ²⁾	+ ²⁾	+ ²⁾	+ ²⁾

Explanation of symbols: + applicable - not applicable ○ limited application

¹⁾ In order to avoid mistakes with regards to service, operation and/or when loading the sterilizer, please refer to table 1 for correct or wrong combination of sterilization containers for standardized steam sterilization procedures.

²⁾ Suitable are only sterilization containers guaranteeing evaporation or discharge of condensate.

³⁾ Gravitation procedures having no drying device do not guarantee a sufficient after-drying effect.

basic model



Basic model containers are produced with perforated or non-perforated bottom. Filters are not protected against contamination from liquids, dust and other mechanical effects. Lids are available in 6 different colours.



General front view



Each container contains two identification labels. One indicating the surgical department where the container belongs; one for the name of the operation. They are easily slid into their built-in slots and firmly secured by a spring enforcement.




Information on the indicator label are filled in by an authorized surgeon before the sterilisation process. These labels are inserted into the spring reinforced slots.




The security seal keeps the container safely locked to avoid any unauthorized opening. The seal automatically breaks when locked flap is lifted.




Safe containers are fitted with a safety lid protecting filters. They are supplied either with filters or without filters. Safety lids are supplied in 6 different colours.



Locking flaps (1) are movable 180° on a pivot and can be secured in any position. Silicon sealings (2) are provided around the lid, in the center and periphery of the round filter retainer, completely recessed.



The filter retainer is locked by pressing down with a thumb in the direction of (1) and open when pressed with a thumb in the direction of (2).




To lift off container lid: push lock (4) in the direction of arrow and hold lid latch, pull in the direction of arrow (5). For mounting safety lid, simply place it on container lid and press down.

1/1 size containers



Types and Dimensions

Bottom Non Perforated



Bottom Perforated



580 x 280 x 100 mm



- 70-7110-10
- 70-7120-10
- 70-7130-10
- 70-7140-10
- 70-7150-10
- 70-7160-10

- 70-7111-10
- 70-7121-10
- 70-7131-10
- 70-7141-10
- 70-7151-10
- 70-7161-10

580 x 280 x 135 mm



- 70-7110-13
- 70-7120-13
- 70-7130-13
- 70-7140-13
- 70-7150-13
- 70-7160-13

- 70-7111-13
- 70-7121-13
- 70-7131-13
- 70-7141-13
- 70-7151-13
- 70-7161-13

580 x 280 x 150 mm



- 70-7110-15
- 70-7120-15
- 70-7130-15
- 70-7140-15
- 70-7150-15
- 70-7160-15

- 70-7111-15
- 70-7121-15
- 70-7131-15
- 70-7141-15
- 70-7151-15
- 70-7161-15

580 x 280 x 200 mm



- 70-7110-20
- 70-7120-20
- 70-7130-20
- 70-7140-20
- 70-7150-20
- 70-7160-20

- 70-7111-20
- 70-7121-20
- 70-7131-20
- 70-7141-20
- 70-7151-20
- 70-7161-20

580 x 280 x 260 mm



- 70-7110-26
- 70-7120-26
- 70-7130-26
- 70-7140-26
- 70-7150-26
- 70-7160-26


- 70-7111-26
- 70-7121-26
- 70-7131-26
- 70-7141-26
- 70-7151-26
- 70-7161-26

1/1 size containers


safe model

Types and Dimensions

Bottom Non Perforated



Bottom Perforated



580 x 280 x 100 mm



- 70-7112-10
- 70-7122-10
- 70-7132-10
- 70-7142-10
- 70-7152-10
- 70-7162-10

- 70-7113-10
- 70-7123-10
- 70-7133-10
- 70-7143-10
- 70-7153-10
- 70-7163-10

580 x 280 x 135 mm



- 70-7112-13
- 70-7122-13
- 70-7132-13
- 70-7142-13
- 70-7152-13
- 70-7162-13

- 70-7113-13
- 70-7123-13
- 70-7133-13
- 70-7143-13
- 70-7153-13
- 70-7163-13

580 x 280 x 150 mm



- 70-7112-15
- 70-7122-15
- 70-7132-15
- 70-7142-15
- 70-7152-15
- 70-7162-15

- 70-7113-15
- 70-7123-15
- 70-7133-15
- 70-7143-15
- 70-7153-15
- 70-7163-15

580 x 280 x 200 mm



- 70-7112-20
- 70-7122-20
- 70-7132-20
- 70-7142-20
- 70-7152-20
- 70-7162-20

- 70-7113-20
- 70-7123-20
- 70-7133-20
- 70-7143-20
- 70-7153-20
- 70-7163-20

580 x 280 x 260 mm



- 70-7112-26
- 70-7122-26
- 70-7132-26
- 70-7142-26
- 70-7152-26
- 70-7162-26

- 70-7113-26
- 70-7123-26
- 70-7133-26
- 70-7143-26
- 70-7153-26
- 70-7163-26

1/1 size containers


bottoms and lids

Types and Dimensions

Bottom Non Perforated



Bottom Perforated



580 x 280 x 100 mm	70-7105-10
580 x 280 x 135 mm	70-7105-13
580 x 280 x 150 mm	70-7105-15
580 x 280 x 200 mm	70-7105-20
580 x 280 x 260 mm	70-7105-26



580 x 280 x 100 mm	70-7106-10
580 x 280 x 135 mm	70-7106-13
580 x 280 x 150 mm	70-7106-15
580 x 280 x 200 mm	70-7106-20
580 x 280 x 260 mm	70-7106-26

Types and Dimensions

Lid Perforated



Safety Lid




580 x 280 mm



70-7100-01
70-7100-02
70-7100-03
70-7100-04
70-7100-05
70-7100-06

580 x 280 mm



70-7103-01
70-7103-02
70-7103-03
70-7103-04
70-7103-05
70-7103-06

accessories

Special Accessories

Wire Baskets



485 x 255 x 30 mm	70-7190-03
485 x 255 x 50 mm	70-7190-05
485 x 255 x 70 mm	70-7190-07
485 x 255 x 100 mm	70-7190-10


Wire Baskets



540 x 255 x 30 mm	70-7191-03
540 x 255 x 50 mm	70-7191-05
540 x 255 x 70 mm	70-7191-07
540 x 255 x 100 mm	70-7191-10

Special Accessories

Wire Base



545 x 255 x 30 mm **70-7193-00**

Drape Retainer




545 x 255 mm **70-7195-00**

Drape



130 x 60 cm	70-7095-07
140 x 100 cm	70-7095-08
185 x 45 cm	70-7095-09
185 x 60 cm	70-7095-10

Silicon Mat



520 x 230 mm	70-7852-30
440 x 230 mm	70-7852-20

3/4 size containers




3/4 size containers


basic model

Types and Dimensions

Bottom Non Perforated



Bottom Perforated



465 x 280 x 100 mm



- 70-7210-10
- 70-7220-10
- 70-7230-10
- 70-7240-10
- 70-7250-10
- 70-7260-10

- 70-7211-10
- 70-7221-10
- 70-7231-10
- 70-7241-10
- 70-7251-10
- 70-7261-10

465 x 280 x 135 mm



- 70-7210-13
- 70-7220-13
- 70-7230-13
- 70-7240-13
- 70-7250-13
- 70-7260-13

- 70-7211-13
- 70-7221-13
- 70-7231-13
- 70-7241-13
- 70-7251-13
- 70-7261-13

465 x 280 x 150 mm



- 70-7210-15
- 70-7220-15
- 70-7230-15
- 70-7240-15
- 70-7250-15
- 70-7260-15


- 70-7211-15
- 70-7221-15
- 70-7231-15
- 70-7241-15
- 70-7251-15
- 70-7261-15

3/4 size containers


safe model

Types and Dimensions

Bottom Non Perforated



Bottom Perforated



465 x 280 x 100 mm



- 70-7212-10
- 70-7222-10
- 70-7232-10
- 70-7242-10
- 70-7252-10
- 70-7262-10

- 70-7213-10
- 70-7223-10
- 70-7233-10
- 70-7243-10
- 70-7253-10
- 70-7263-10

465 x 280 x 135 mm



- 70-7212-13
- 70-7222-13
- 70-7232-13
- 70-7242-13
- 70-7252-13
- 70-7262-13

- 70-7213-13
- 70-7223-13
- 70-7233-13
- 70-7243-13
- 70-7253-13
- 70-7263-13

465 x 280 x 150 mm



- 70-7212-15
- 70-7222-15
- 70-7232-15
- 70-7242-15
- 70-7252-15
- 70-7262-15

- 70-7213-15
- 70-7223-15
- 70-7233-15
- 70-7243-15
- 70-7253-15
- 70-7263-15

3/4 size containers

bottoms and lids

Types and Dimensions

Bottom Non Perforated



Bottom Perforated



465 x 280 x 100 mm **70-7205-10**
465 x 280 x 135 mm **70-7205-13**
465 x 280 x 150 mm **70-7205-15**



465 x 280 x 100 mm **70-7206-10**
465 x 280 x 135 mm **70-7206-13**
465 x 280 x 150 mm **70-7206-15**

Types and Dimensions

Lid Perforated



Safety Lid



70-7200-01
70-7200-02
70-7200-03
70-7200-04
70-7200-05
70-7200-06



70-7203-01
70-7203-02
70-7203-03
70-7203-04
70-7203-05
70-7203-06

3/4 size containers accessories


Special Accessories

Wire Baskets



405 x 255 x 30 mm **70-7290-03**
405 x 255 x 50 mm **70-7290-05**
405 x 255 x 70 mm **70-7290-07**
405 x 255 x 100 mm **70-7290-10**

Wire Base



405 x 255 x 30 mm **70-7293-03**


Special Accessories

Drape Retainer




405 x 255 mm **70-7295-00**

Drape



140 x 100 cm **70-7095-08**

Silicon Mat



440 x 230 mm **70-7852-20**
380 x 230 mm **70-7852-10**

1/2 size containers




1/2 size containers

basic model

Types and Dimensions

Bottom Non Perforated



Bottom Perforated



285 x 280 x 100 mm



- 70-7310-10
- 70-7320-10
- 70-7330-10
- 70-7340-10
- 70-7350-10
- 70-7360-10

- 70-7311-10
- 70-7321-10
- 70-7331-10
- 70-7341-10
- 70-7351-10
- 70-7361-10

285 x 280 x 135 mm



- 70-7310-13
- 70-7320-13
- 70-7330-13
- 70-7340-13
- 70-7350-13
- 70-7360-13

- 70-7311-13
- 70-7321-13
- 70-7331-13
- 70-7341-13
- 70-7351-13
- 70-7361-13

285 x 280 x 150 mm



- 70-7310-15
- 70-7320-15
- 70-7330-15
- 70-7340-15
- 70-7350-15
- 70-7360-15

- 70-7311-15
- 70-7321-15
- 70-7331-15
- 70-7341-15
- 70-7351-15
- 70-7361-15

285 x 280 x 200 mm



- 70-7310-20
- 70-7320-20
- 70-7330-20
- 70-7340-20
- 70-7350-20
- 70-7360-20

- 70-7311-20
- 70-7321-20
- 70-7331-20
- 70-7341-20
- 70-7351-20
- 70-7361-20

285 x 280 x 260 mm




- 70-7310-26
- 70-7320-26
- 70-7330-26
- 70-7340-26
- 70-7350-26
- 70-7360-26


- 70-7311-26
- 70-7321-26
- 70-7331-26
- 70-7341-26
- 70-7351-26
- 70-7361-26

Types and Dimensions

Bottom Non Perforated



Bottom Perforated



285 x 280 x 100 mm



-  70-7312-10
-  70-7322-10
-  70-7332-10
-  70-7342-10
-  70-7352-10
-  70-7362-10

-  70-7313-10
-  70-7323-10
-  70-7333-10
-  70-7343-10
-  70-7353-10
-  70-7363-10

285 x 280 x 135 mm




-  70-7312-13
-  70-7322-13
-  70-7332-13
-  70-7342-13
-  70-7352-13
-  70-7362-13

-  70-7313-13
-  70-7323-13
-  70-7333-13
-  70-7343-13
-  70-7353-13
-  70-7363-13

285 x 280 x 150 mm



-  70-7312-15
-  70-7322-15
-  70-7332-15
-  70-7342-15
-  70-7352-15
-  70-7362-15

-  70-7313-15
-  70-7323-15
-  70-7333-15
-  70-7343-15
-  70-7353-15
-  70-7363-15

285 x 280 x 200 mm



-  70-7312-20
-  70-7322-20
-  70-7332-20
-  70-7342-20
-  70-7352-20
-  70-7362-20

-  70-7313-20
-  70-7323-20
-  70-7333-20
-  70-7343-20
-  70-7353-20
-  70-7363-20

285 x 280 x 260 mm



-  70-7312-26
-  70-7322-26
-  70-7332-26
-  70-7342-26
-  70-7352-26
-  70-7362-26

-  70-7313-26
-  70-7323-26
-  70-7333-26
-  70-7343-26
-  70-7353-26
-  70-7363-26

Bottom and Lids


Types and Dimensions

Bottom Non Perforated



285 x 280 x 100 mm	70-7305-10
285 x 280 x 135 mm	70-7305-13
285 x 280 x 150 mm	70-7305-15
285 x 280 x 200 mm	70-7305-20
285 x 280 x 260 mm	70-7305-26

Bottom Perforated



285 x 280 x 100 mm	70-7306-10
285 x 280 x 135 mm	70-7306-13
285 x 280 x 150 mm	70-7306-15
285 x 280 x 200 mm	70-7306-20
285 x 280 x 260 mm	70-7306-26

Types and Dimensions

Lid Perforated



285x 280 mm



70-7300-01
70-7300-02
70-7300-03
70-7300-04
70-7300-05
70-7300-06

Safety Lid



285 x 280 mm



70-7303-01
70-7303-02
70-7303-03
70-7303-04
70-7303-05
70-7303-06

1/2 size containers accessories


Special Accessories

Wire Baskets



255 x 245 x 30 mm	70-7390-03
255 x 245 x 50 mm	70-7390-05
255 x 245 x 70 mm	70-7390-07
255 x 245 x 100 mm	70-7390-10

Wire Base



258 x 258 x 30 mm	70-7393-00
-------------------	-------------------

Special Accessories

Drape Retainer



255 x 255 mm	70-7395-00
--------------	-------------------

Drape



90 x 60 cm	70-7095-04
90 x 80 cm	70-7095-05
130 x 45 cm	70-7095-06
130 x 60 cm	70-7095-17

Silicon Mat



245 x 245 mm	70-7851-00
220 x 230 mm	70-7852-00

other containers




Wide body containers

basic model

Types and Dimensions

Bottom Non Perforated



Bottom Perforated



600 x 400 x 120 mm



-  70-7110-62
-  70-7120-62
-  70-7130-62
-  70-7140-62
-  70-7150-62
-  70-7160-62

-  70-7111-62
-  70-7121-62
-  70-7131-62
-  70-7141-62
-  70-7151-62
-  70-7161-62

600 x 400 x 180 mm



-  70-7110-68
-  70-7120-68
-  70-7130-68
-  70-7140-68
-  70-7150-68
-  70-7160-68

-  70-7111-68
-  70-7121-68
-  70-7131-68
-  70-7141-68
-  70-7151-68
-  70-7161-68

Wide body containers

bottoms, lids and accessories

Types and Dimensions

Bottom Non Perforated



Bottom Perforated



Lid Perforated



70-7100-51

70-7100-52

70-7100-53

70-7100-54

70-7100-55

70-7100-56

600 x 400 x 120 mm
600 x 400 x 180 mm

70-7105-62
70-7105-68

600 x 400 x 120 mm
600 x 400 x 180 mm

70-7106-62
70-7106-88


Special Accessories

Wire Baskets



470 x 337 x 51 mm

70-7190-55




533 x 337 x 101 mm

70-7191-60

Types and Dimensions

Bottom Non Perforated






285 x 280 x 55 mm




285x 280 x 85 mm



-  70-7310-05
-  70-7320-05
-  70-7330-05
-  70-7340-05
-  70-7350-05

-  70-7310-08
-  70-7320-08
-  70-7330-08
-  70-7340-08
-  70-7350-08

Bottom Perforated



-  70-7311-08
-  70-7321-08
-  70-7331-08
-  70-7341-08
-  70-7351-08


Bottoms, lids and accessories

Types and Dimensions

Bottom Non Perforated



Bottom Perforated



Lid Perforated



285 x 280 x 55 mm
285 x 280 x 85 mm

70-7305-05
70-7305-08



285 x 280 x 85 mm **70-7306-08**



70-7304-01

70-7304-02


70-7304-03

70-7304-04

70-7304-05

Special Accessories

Wire Baskets



255 x 245 x 30 mm
255 x 245 x 50 mm

70-7390-03
70-7390-05

Wire Base



245 x 255 x 30 mm **70-7393-00**

Silicon Mat



245 x 245 mm
220 x 230 mm

70-7851-00
70-7852-00

Drape Retainer



255 x 255 mm **70-7395-00**

Drape



90 x 60 cm
90 x 80 cm

70-7095-04
70-7095-05

Types and Dimensions

Bottom Non Perforated



Bottom Perforated



300 x 140 x 40 mm






-  70-7410-04
-  70-7420-04
-  70-7430-04
-  70-7440-04
-  70-7450-04

-  70-7411-04
-  70-7421-04
-  70-7431-04
-  70-7441-04
-  70-7451-04

300 x 140 x 70 mm



-  70-7410-07
-  70-7420-07
-  70-7430-07
-  70-7440-07
-  70-7450-07

-  70-7411-07
-  70-7421-07
-  70-7431-07
-  70-7441-07
-  70-7451-07

300 x 140 x 100 mm



-  70-7410-10
-  70-7420-10
-  70-7430-10
-  70-7440-10
-  70-7450-10

-  70-7411-10
-  70-7421-10
-  70-7431-10
-  70-7441-10
-  70-7451-10

bottoms, lids and accessories

Types and Dimensions

Bottom Non Perforated



Bottom Perforated



Lid Perforated




300 x 140 x 40 mm **70-7405-04**
300 x 140 x 70 mm **70-7405-07**
300 x 140 x 100 mm **70-7405-10**



300 x 140 x 40 mm **70-7406-04**
300 x 140 x 70 mm **70-7406-07**
300 x 140 x 100 mm **70-7406-10**

300 x 140 mm



70-7400-01
70-7400-02
70-7400-03
70-7400-04
70-7400-05

Special Accessories

Silicon Mat



270 x 125 mm **70-7850-00**

Drape Retainer



280 x 125 mm **70-7495-00**

Drape




60 x 40 cm **70-7095-01**

Types and Dimensions

Lid and Bottom Perforated



Bottom Perforated



Lid Perforated



310 x 190 x 40 mm



70-7510-04

310 x 190 x 65 mm




70-7510-06

310 x 190 x 130 mm



70-7510-13



310 x 190 x 40 mm **70-7506-04** 310 x 190 mm **70-7500-01**
 310 x 190 x 65 mm **70-7506-06**
 310 x 190 x 130 mm **70-7506-13**

Special Accessories

Drape Retainer



310 x 190 mm **70-7595-00**

Drape



60 x 40 cm **70-7095-01**
 65 x 65 cm **70-7095-02**
 85 x 85 cm **70-7095-03**

Types and Dimensions

Lid Perforated - Bottom Perforated (Non Filter)



300 x 138 x 65 mm



- 70-7610-06
- 70-7620-06
- 70-7630-06
- 70-7640-06
- 70-7650-06

500 x 155 x 75 mm



- 70-7612-07
- 70-7622-07
- 70-7632-07
- 70-7642-07
- 70-7652-07

Scopy containers

Types and Dimensions

Lid Perforated - Bottom Perforated (Non Filter)



250 x 60 x 50 mm



- 70-7600-01

450 x 70 x 70 mm




- 70-7600-02

Removal containers

Bottom, lids and accessories

Types and Dimensions

Lid And Bottom Non Perforated



Bottom Non Perforated



285 x 280 x 100 mm



285 x 280 x 150 mm



580 x 280 x 150 mm



- 70-7315-10
- 70-7325-10
- 70-7335-10
- 70-7345-10
- 70-7355-10
- 70-7365-10

- 70-7315-15
- 70-7325-15
- 70-7335-15
- 70-7345-15
- 70-7355-15
- 70-7365-15

- 70-7115-15
- 70-7125-15
- 70-7135-15
- 70-7145-15
- 70-7155-15
- 70-7165-15



285 x 280 x 100 mm **70-7305-10**

285 x 280 x 150 mm **70-7305-15**



580 x 280 x 150 mm **70-7105-15**

Lids - Non Perforated



Lids - Non Perforated



580 x 280 mm




- 70-7102-01
- 70-7102-02
- 70-7102-03
- 70-7102-04
- 70-7102-05
- 70-7102-06

285 x 280 mm



- 70-7302-01
- 70-7302-02
- 70-7302-03
- 70-7302-04
- 70-7302-05
- 70-7302-06



70-7880-00

Brush container

285 x 280 x 85 mm

Containers are supplied together with mounting fixtures onto the wall.

For 20 brushes. Brushes are taken out by moving of the handle in both directions right or left which is under the container.



70-7880-00

70-7880-01

70-7880-02

70-7880-03

70-7880-04

Filters, Indicators And Labels


Disposable paper filters



For 1/1, 3/4 and Half Size Containers ø 190 mm
70-7800-00
 100 Pcs/Package
70-7800-00G
 1000 Pcs/Package



70-7800-10
 100 Pcs/Package
70-7800-10G
 1000 Pcs/Package



70-7802-00
 100 Pcs/Package
 For Mini Containers

70-7804-00
 100 Pcs/Package
 For Dental Containers




70-7801-00
 Filter Dispenser

Textile Filters



70-7800-01
 Single
 For 1/1, 3/4 and Half Size Containers



70-7802-01
 Single
 For Mini Containers
70-7804-01
 Single
 For Dental Containers



70-7809-00
 For Paper Filter
70-7809-01
 For Textile Filter




70-7810-00
 100 Pcs / Package
 Paper label with indicator
 For 1/1, 3/4 and Half Size Containers



70-7812-00
 100 Pcs / Package
 Paper label with indicator
 For Mini and Dental Containers




70-7819-00
 100 Pcs / Package
 Plastic Security Seal



- 70-7820-01
- 70-7820-02
- 70-7820-03
- 70-7820-04
- 70-7820-05
- 70-7820-06

Coloured Identification Labels
 Engraved to your specifications, max. 16 letters.



70-7829-00
 Dispenser for identification labels

Components



70-7830-00
Holding Pin
For Wire Baskets




70-7830-01
Holding Pin
For Partition Sheets



Partition Sheets
Stainless Steel


- 70-7831-05** 50 x 20 mm
- 70-7831-13** 130 x 20 mm
- 70-7831-22** 225 x 20 mm
- 70-7831-46** 460 x 20 mm




70-7840-00
Divider, to be used with 70-7841-00



70-7841-00 232 mm
70-7841-01 471 mm
Sliding Rails with Securing Screws




70-7842-00 232 mm
70-7842-01 471 mm
Limiting Bar with Securing Screws



70-7845-00
Spiral Holding



70-7847-00
Double Spiral Holding



Fixation clamps



70-7849-00 50 x 20 x 12
70-7849-01 65 x 20 x 20



70-7849-20 Ø 7-12 mm
70-7849-21 Ø 12-16 mm
70-7849-22 Ø 16-26 mm
70-7849-23 Ø 26-36 mm
70-7849-24 Ø 36-45 mm




Carrying



70-7890-00
Trolley
720 x 620 x 800 mm



70-7891-00
Trolley
1030 x 620 x 1170 mm



Distributed by:

Art-N° 90-1215



 **SCHREIBER**[®] INSTRUMENTE

SCHREIBER GmbH
Unterer Damm 15
D-78567 Fridingen, Germany
Tel.: +49 7463 996949-0
Fax.: +49 7463 996949-9
E-Mail: info@schreibergmbh.de
Internet: www.schreibergmbh.de