



ULTRABUR OTOLOGIC DRILL

- Hand-held microdrill for stapedotomy, stapedectomy, and fashion of homo/auto bone grafts
- Narrow, lightweight design allows for better visualization in middle ear
- Speed: 7,000 rpm max
- Provided sterile

Number			Qty/Box
10-3030	Diamond tip	0.6 mm	1
10-3031	Diamond tip	0.7 mm	1
10-3038	Diamond tip	0.8 mm	1

BURS

- Designed for use with straight shank drill systems and Xomed® power hand pieces
- Multiple uses
- Provided sterile

Round Carbide

Number w/ Notch	w/o Notch	Head Diameter (mm)	Length (mm)	Number of Flutes
BC-70050		0.5	70	4
BC-70100	BC-70100P	1.0	70	6
BC-70150		1.5	70	6
BC-70200	BC-70200P	2.0	70	6
BC-70250		2.5	70	6
BC-70300	BC-70300P	3.0	70	8
BC-70350		3.5	70	8
BC-70400	BC-70400P	4.0	70	8
BC-70450		4.5	70	10
BC-70500	BC-70500P	5.0	70	10
BC-70550		5.5	70	10
BC-70600	BC-70600P	6.0	70	10
BC-70650		6.5	70	10
BC-70700	BC-70700P	7.0	70	12
BC-70000	BC-70000P	Complete Set - one of each		

Round Diamond Bur

Number w/ Notch	w/o Notch	Head Diameter (mm)	Length (mm)	Style
BD-70050		0.5	70	Diamond
BD-70100	BD-70100	1.0	70	Diamond
BD-70150		1.5	70	Diamond
BD-70200	BD-70200P	2.0	70	Diamond
BD-70250		2.5	70	Diamond
BD-70300	BD-70300P	3.0	70	Diamond
BD-70350		3.5	70	Diamond
BD-70400	BD-70400P	4.0	70	Diamond
BD-70500	BD-70500P	5.0	70	Diamond
BD-70600		6.0	70	Diamond
BD-70000	BD-70000P	Complete Set - one of each		
O225300		Bur Rack - Holds 34 straight shaft burs, length 70 mm		