


# Surg-I-Loop®

## Silicone Loops



Length  
16 1/2" (42 cm)

Designed to assist surgeons in complicated surgical fields by providing retraction, occlusion and identification of arteries, veins, nerves and ureters.

- **Unique design** – provides access to complicated surgical fields
- **Versatile** – can retract and occlude arteries, veins, nerves and ureters
- **Durable** – retains preset tension
- **Radiopaque**
- **Latex-Safe**

- Super Maxi (5 x 1.27 mm)
- Maxi (2.5 x 1.02 mm)
- Mini (1.5 x 0.88 mm)
- Micro (0.5 x 0.26 mm)

Cat. No.	Description
1001-75	Red, Mini
1001-76	Red, Maxi
1001-77	Blue, Mini
1001-78	Blue, Maxi
1001-78A	Blue, Super Maxi
1001-79	Yellow, Mini
1001-80	Yellow, Maxi
1001-81	White, Mini
1001-81M	White, Micro
1001-82	White, Maxi

*Two per sterile package,  
10 packages per box*