



Designed for the surgical management of tracheal and subglottic stenosis and the reconstruction of cervical and thoracic trachea. Packaged 1 each.

- 100% medical grade silicone
- Easily trimmed to meet patient's needs
- Large soft collar with closure attached

Catalog No.	Size	A	B	C	D	E	F
70-50080	Standard	8mm	37mm	16.5mm	32mm	6.5mm	14mm
70-50100	Standard	10mm	50mm	24mm	42mm	8mm	18mm
70-50120	Standard	12 m	66mm	32mm	50mm	10.5mm	23mm
70-50140	Standard	14mm	71mm	34mm	50mm	10.5mm	26mm
70-50160	Standard	16mm	72mm	34mm	60mm	11.5mm	26mm
70-50081	Long	8mm	80mm	60mm	32mm	6.5mm	14mm
70-50101	Long	10mm	96mm	70mm	42mm	8mm	18mm
70-50121	Long	12mm	108mm	75mm	50mm	10.5mm	23mm
70-50141	Long	14mm	112mm	75mm	50mm	10.5mm	26mm
70-50161	Long	16mm	122mm	85mm	60mm	11.5mm	26mm