

TUBOS DE PENROSE

"LÁTEX - NULIFE[®] - X RAY"



Dimensiones y medidas:

REFERENCIA	DESCRIPCIÓN	ANCHO TUBO	Ø INTERIOR	LONGITUD	CAJA	EMBALAJE
27.04.09	TUBO PENROSE NF LATEX	9.5mm	Ø 1/4"	±45 cm.	50	800
27.04.13	TUBO PENROSE NF LATEX	13mm	Ø 3/8"	±45 cm.	50	800
27.04.19	TUBO PENROSE NF LATEX	19mm	Ø 1/2"	±45 cm.	50	800
27.04.25	TUBO PENROSE NF LATEX	25mm	Ø 5/8"	±45 cm.	50	800
27.04.32	TUBO PENROSE NF LATEX	32mm	Ø 7/8"	±45 cm.	50	800
27.04.38	TUBO PENROSE NF LATEX	38mm	Ø 1"	±45 cm.	50	800

1 pulgada = 25'4mm

Presentación: Envasados individualmente y estériles.